

# EnFix<sup>®</sup>

DEMINERALIZED BONE FIBER IMPLANT FOR ENTESIS REPAIR

EnFix RC<sup>™</sup>



EnFix TAC<sup>®</sup>

IT'S ALL ABOUT THE ENTESIS<sup>®</sup>

# EnFix<sup>®</sup>

DEMINERALIZED BONE FIBER IMPLANT FOR ENTHESIS REPAIR

## Enthesis Failure Syndrome


Lack of healing at the tendon bone interface is a common mode of failure in rotator cuff repair.

### Rotator Cuff Repair

In the United States alone, nearly 500,000 rotator cuff repairs are performed annually and 20% to 70% of these repairs fail structurally. Inadequate tendon-to-bone ingrowth results in incomplete healing, gap formation and a higher risk of re-tear. Tendon reattachment is a crucial clinical need, especially in larger tears because failure rates increase linearly with tear size.

Most augmentation products have been designed as “overlays” to reinforce the tendon. EnFix RC and EnFix TAC change the paradigm by enhancing healing at the enthesis where failure often occurs. This enhanced biologic repair at the interface from the bone to the tendon is a significant advance. EnFix implants are produced using demineralized bone fiber (DBF) to provide optimal biologic performance while also easily integrating into current surgical techniques.

The images below demonstrate reformation of the enthesis at 12 weeks following treatment with EnFix RC in a pre-clinical sheep study.



H&E histology at 12 weeks under normal light (top left and right) for DBF treated enthesis revealed an active interface with some residual DBF and enthesis reformation. Polarized light (bottom) confirmed Sharpey's fibers in the DBF treated groups that were not present in the controls.


## Surgical Technique

### Preparing the Implant Site

- The EnFix products interface easily with existing surgical techniques.
- Prepare implant site as usual.
- Use an awl to make the insertion site.


### EnFix RC Surgical Technique

- Place device on introducer and insert the graft fully so the top of the device is flush with, or on the cortical surface.
- Tap device to receive a suture anchor using suture anchor's tap.
- Insert suture anchor.
- Insert additional devices / suture anchors as required.
- Complete repair in usual manner.



### EnFix TAC Surgical Technique

- Place device on introducer and insert the graft fully so the top of the device is flush with, or on the cortical surface.
- Insert additional devices as required.
- Complete repair in usual manner.




# EnFix®

DEMINERALIZED BONE FIBER IMPLANT FOR ENTESIS REPAIR

## Rabbit Infrapinatus Repair

In a rabbit infrapinatus model, the tendon was detached and reattached using suture fixation through two bone tunnels. Demineralized bone fiber was placed between the tendon and the bone in the treated group. Sutures were cut prior to mechanical testing at each time interval.


The graph below shows increased strength with the treated group vs. the control, and the control was not statistically higher at 12 weeks vs. 6 weeks.



## Suture Anchor Pull Out Testing

Pull out strength testing was performed in Sawbones Foam (1522-09; 10 pfc), a bone analog specified in ASTM standards for pull out testing. In the treated group, an EnFix RC implant (4.5mm) was inserted in the usual manner and a 5.5mm PEEK suture anchor was gently tapped into the device. The control group used a suture anchor alone. Pull out testing was performed at 20 mm/min, with maximum load recorded.


The graph below shows increased pull out force with the EnFix RC implant vs. the control.



## FormLok™ Technology

The EnFix family of implants are 100% cortical bone and there are no excipients, therefore they conform with the regulations for minimally manipulated tissue. The FormLok process imparts shape retention to the device, even when immersed in liquid, as is often required for use in arthroscopic surgery.


The figures below show timelapse images of EnFix with and without FormLok treatment. The nontreated sample on the left rapidly loses its shape, while the FormLok treated device on the right retains its shape at 15 minutes, and beyond.




## Nanotopography

Topography of a surface can influence cellular response. The process to manufacture the DBF fibers is designed to conserve the collagen structure of bone while preserving the inductive proteins (the BMPs). This is achieved by the proprietary process that demineralizes the bone and the fibers are then created by cleaving the bone along the axis of collagen orientation along the surfaces of the collagen fibrils. This gentle process provides a nanotopography that is not seen in conventional bone matrix products that are acid treated to demineralize the bone after particle or fiber formation.

The DBF fiber nanotopography of the Tetrous fibers is shown in the top images below. In contrast, the lower images shown below are of a conventional DBM particle where it can be seen that the acid treatment smooths the surface destroying its beneficial topography.



Tetrous DBF fibers at low and high magnification



Conventional DBM at the same magnifications

# EnFix®


DEMINERALIZED BONE FIBER IMPLANT FOR ENTESIS REPAIR



EnFix RC



EnFix TAC-T



EnFix TAC-O

<u>Part Number</u>	<u>Size/Description</u>
TET-RC-45	For use with suture anchors 4.5mm to 5.5mm
TET-RC-55	For use with suture anchors larger than 5.5mm

<u>Part Number</u>	<u>Size/Description</u>
TET-TAC-T	4mm x 10mm

For use with all-suture anchor or independent of an anchor

<u>Part Number</u>	<u>Size/Description</u>
TET-TAC-O	8.5mm diameter

TETROUS, INC.  
2021 MIDWEST RD., STE. 101  
OAK BROOK, IL 60523

1-331-307-7499  
[www.tetrous.com](http://www.tetrous.com)

Tetrous®, EnFix®, EnFix RC™, EnFix TAC® and It's All About The Entesis® are trademarks of Tetrous, Inc. FormLok™ is a trademark of TheraCell, Inc., an ISTO Biologics company.

US 9,486,557, US 9,572,912, US 11,660,373,  
US 11,759,548, US 12,036,338, US 12,115,279,  
US 12,171,908  
Other patents pending.  
For more information, see [www.tetrous.com/patents](http://www.tetrous.com/patents).

The Tetrous DBF technology and other processes are licensed exclusively from TheraCell, Inc., an ISTO Biologics company.

ALL RIGHTS RESERVED

© 2025 Tetrous, Inc.

<u>Part Number</u>	<u>Description</u>	<u>Size</u>
MO-AWL-0650-A	Standard Awl	4.5mm/5.5mm
MO-AWL-0650-C	Cutting Awl	4.5mm
MO-AWL-0650-B	Cutting Awl	5.5mm
MO-PRB-0651-A	Introducer	
MO-TRA-0670	Instrument Tray	4 Instruments

